# Nite Tibial Nail

Nite is a titanium nail designed for the treatment of tibial fractures. By using Claw technology for distal locking, Nite helps healthcare professionals to achieve superior operative and postoperative parameters.

Nite has a fully cannulated body compatible with guide-wire applications and offers 7 mm compression range, giving the surgeons more resources during the operation.

#### Indications

- Proximal extra-articular fracture
- Open and closed fractures of the tibial shaft
- Pathologic / impending pathologic fractures
- Malunions / nonunions

## Nail's key figures

#### Nail length

270 mm to 375 mm in 15 mm increment

- Proximal Diameter
  11 mm
- Distal Diameter

9 mm to 11 mm in 1 mm increment

- Distal Claw Maximum Span
  - 38 mm
- Compression Range7 mm
- End Cap Length
  - 0 mm to 10 mm in 5 mm increments Internal thread to secure the cap to the 5mm Hex driver

### Compression and Cortical Screws' key figures

Length

 $\checkmark$ 

30 mm to 120 mm in 5 mm increment

- Diameter
  - 5 mm

Internal thread to secure the screw to the 5mm Hex driver. Compression Screw with threaded tip and 5 mm shaft to withstand the compression load Claws are deployed from within the nail, and do not protrude from the bone. This eliminates the **risk of vascular or nerve damage** as well as **discomfort** caused by conventional screws on the distal tibia.